


Curatest®

Epikutantestpflaster
Emplâtre de test épicutané
Patch test strip
Apósito para pruebas epicutáneas
Adesivo para testes epicutâneos
Cerotto per test epicutanei
Epicutaantestpleister
Epikutantestplaster
Epikutantestplåster
Epikutaanitestilaastari
Náplast pro epikutánní testy
Epikutánná testovacia náplast
Epikután próba tapasz
Plaster do testów skórnych
Oblič za epikutano testiranje
Flaster za epikutani test

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MD Dispositif médical

en Instructions for use

Product description and features:
White strip that can be applied with test substances to perform epicutaneous tests. The strip consists of an adhesive-coated nonwoven material. It comprises a separately detachable middle section with 10 integrated test discs and two nonwoven strips (margins).

Product composition:
Polypropylene, polyacrylate, viscose, polyethylene, polyester, silicone paper

Intended purpose:
The Curatest® Adhesive Strip for Patch Test is intended to be used for patch-tests.

Indications:
Not applicable.

Contraindications:
Known allergy and/or hypersensitivity to any of the product components.

Side effects:
In very rare cases skin irritations and/or allergies may occur.

Application instructions:
The product is applied by healthcare professionals.

- Preparation:
Before application, please ensure that the test area is free of hair and that the skin is dry and grease-free to ensure that the patch test strip adheres properly and that there is contact between the skin and the allergens.
- Use:
 - Remove a Curatest® patch test strip from the folding carton and label the margins on the nonwoven side (Fig. 1).
 - Open the two cover papers up to the perforation (Fig. 2).
 - Coat the test discs with the test substances (Fig. 3).
 - Take hold of the opened patch test strip with your hand so as to not touch the adhesive surface of the strip and apply with the adhesive side on the desired part of the skin (Fig. 4).
 - Open the two cover papers completely. For secure fixation, pass your hand particularly over the edges of the patch test strip (Fig. 5).
 - Follow-up appointment: Remove the middle section of the patch test strip and leave the two margins in place. Possible skin irritations can thus be matched with the labeled test substances (Fig. 6).
- Disposal:
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions:
In case of a serious incident, contact the manufacturer and the responsible health authorities.

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