Rosidal® K Rosidal® SC Elko® - Rosidal®

50x

Durelast® Durelast® Combi

10x

Rosidal® haft

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en Instructions for Use

Product description and performance features:

- Rosidal® K / Elko® Rosidal®: stretch approx. 90%
- Rosidal® SC: stretch approx. 50%
- Durelast® / Durelast® Combi: stretch approx. 40%
- Rosidal® haft: stretch approx. 60%
- Rosidal® CC: stretch approx. 75%

Product composition:

- Rosidal® K / Elko® Rosidal®: 100% cotton
- Rosidal® SC: 60% polyamide, 40% cotton
- Durelast® / Durelast® Combi: 66% cotton, 34% polyamide
- Rosidal® haft: 100% cotton, with a natural rubber latex microcoating
- Rosidal® CC: 51% cotton, 39% polyamide, 10% polyurethane, with a synthetic microcoating

Intended use:

For compression, retention of dressings, support and relief and for immobilisation of limbs.

Indications:

phlebological and lymphological indications: varicosis, venous insufficiency, thromboangiitis, phlebothrombosis, post-thrombotic syndrome, stage I to III b chronic venous insufficiency according to Widmer, venous and mixed ulcers, edema (phleboedema, lymphoedema,

lypoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery, contusions, sprains, tendon injuries and joint dislocations

Contraindications:

Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects:

Skin irritations and/or allergies can occur in very rare cases.

Warnings and precautions:

Incorrect bandaging can result in skin necrosis and pressure damage to nerves. Slight bluish tinge of the toes is normal, this should disappear as soon as the patient starts to walk. If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression. The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage. For use with ABPIs over 0.8. When ABPI is between 0.5 and 0.8, apply under medical supervision only. Use adhesive strips to fix the bandage in place. Do not fix in place with bandage clips because these can cause injury.

General instructions:

- Rosidal® SC – apply the soft compression bandage with a 50% overlap (C1-C3).
- Rosidal® CC – apply the bandage at full stretch (Safe-Loc®) in figures of 8 (images D1-D3).
- The bandaging technique is chosen in consultation with the treating physician.
- It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.
- The width of the bandage should correspond to the diameter of the limb to be bandaged.
- To achieve maximum joint movement, it is recommended to apply the bandage with the foot at an angle of 90° to the leg.
- The bandage must be applied so that the pressure from the bandage decreases from the foot to the knee and, where applicable, to the gluteal fold.
- The bandage is unrolled with the rolled-up portion on top. Hold the top of the bandage close to the leg/arm.

Usage information:
• Rosidal® K / Elko® - Rosidal® / Durelast® can be applied using different compression techniques (images A1-A11 or B1-B10).
• Unroll approx. 20 cm from the Rosidal® haft roll prior to application.

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