Rosidal® Lymph

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en Instructions for Use

Product description and performance characteristics

Compression system for the treatment of lymphological indications.

Product composition

- Mollelast®: 56 % viscose, 44 % polyamide
- Rosidal® SC: polyurethane foam with woven fabric made from 60 % polyamide, 40 % cotton
- Komprex®: 80 % styrene-butadiene rubber
- Rosidal® K: 100 % cotton
- Idealbinde®: 100 % cotton

The number of components varies in the individual Rosidal® Lymph versions (see table).

Intended purpose

For compression

Indications

Primary or secondary lymphoedema, lipoeedema, oedema of the upper and/or lower extremities

Contraindications

Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects

In rare cases skin irritations and/or allergies may occur.

Warnings and precautions

1. Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
2. Slight bluish tinge to the toes is normal. This should disappear as soon as the patient starts to walk.
3. If the symptoms (pain and bluish tinge) do not disappear, the bandage must be removed and applied again with less compression.
4. The patient must inform their doctor

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if sharply increasing pain occurs while wearing the bandage.

2. Then apply the Rosidal® SC soft compression bandage with textile bottom side for gentle compression and to protect the skin and minimise slippage.
3. For even pressure distribution, sections can be cut to fit from the Komprex® foam rubber bandage – suitable for bony protuberances or concave spots – and integrated into the bandage. The edges should be slightly angled to avoid pressure points.
4. Apply the Rosidal® K short-stretch compression bandage on over this. Its effect is due to the working pressure exerted by movement. The patient should remain active despite the compression bandage. For compression treatment of larger areas, the elastic textile Idealbinde® can also be used (only available in Rosidal® Lymph for leg, large).

Recommendations for care

Rosidal® K, Rosidal® SC and Idealbinde®: Recommended washing temperature 40°C–60 °C, can be washed up to 95°C (commercially available detergent), Rosidal® K, Rosidal® SC can be washed up to 50 times, Idealbinde® 5 times. Dry flat.

The maximum stretch of the bandage may increase with washing. Bear this in mind when applying the bandage.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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without tension for gentle compression of the fingers and/or toes.

2. Then apply the Rosidal® SC soft compression bandage with textile bottom side for gentle compression and to protect the skin and minimise slippage.
3. For even pressure distribution, sections can be cut to fit from the Komprex® foam rubber bandage – suitable for bony protuberances or concave spots – and integrated into the bandage. The edges should be slightly angled to avoid pressure points.
4. Apply the Rosidal® K short-stretch compression bandage on over this. Its effect is due to the working pressure exerted by movement. The patient should remain active despite the compression bandage. For compression treatment of larger areas, the elastic textile Idealbinde® can also be used (only available in Rosidal® Lymph for leg, large).

Recommendations for care

Rosidal® K, Rosidal® SC and Idealbinde®: Recommended washing temperature 40°C–60 °C, can be washed up to 95°C (commercially available detergent), Rosidal® K, Rosidal® SC can be washed up to 50 times, Idealbinde® 5 times. Dry flat.

The maximum stretch of the bandage may increase with washing. Bear this in mind when applying the bandage.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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