




Curafix® i.v. Soft

Kanülenfixierpflaster aus Vliesstoff
Pansement de fixation de canule en non-tissé
Cannula Retention Dressing made of nonwoven material
Apósito de fijación de cánulas de tejido sin tejer
Adesivo para fixação de cânulas de tecido não tecido
Cerotto in TnT per il fissaggio delle cannule
Pleister van vliesstof voor het fixeren van canules
Kanylefikseringsplaster af fleece
Nonwoven-kanylfixeringsplåster
Náplast k fixaci kanyl z netkané textilie
Fixačná náplast pre kanylu z netkanej textílie
Kanülrögzítő tapasz flízanyagból
Plaster z włókniny do mocowania kanułu
Obliž iz koprene za pritrdjevanje kanile
Flaster od flizelina za fiksiranje kanila
Dokunmamış kumaştan üretilmiş kanül sabitleme bandı

MD

Medizinprodukt · Dispositif médical · Medical Device · Dispositivo médico · Dispositivo médico · Dispositivo medico · Medisch hulpmiddel · Medientechnisk produkt · Medientechnisk produkt · Zdravotnický prostředek · Zdravotnický prostředek · Zdravotnickai eszköz · Wyrób medyczny · Medicinski pripomoček · Medicinski proizvod · Tibbi ürün

Einfaches Sterilbarriersystem · stème de barrière stérile simple · Single sterile barrier system · Sistema de barrera estéril simple · Sistema de barreira estéril simples · Sistema a singola barriera sterile · Enkelvoudig steriel barrièresysteem · Enkelt sterilt barriersystem · Enkelt sterilbarriärsystem · Jednoduchý sterilní bariérový systém · Egyes steril gátrendszér · Pojedynczy system bariery sterylnej · Enojni sistem sterilne pregrade · Sustav jednostruke sterilne barijere · Tek bariyerli steril bariyer sistemi

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en Instructions for Use

Product description and features:
Curafix® i.v. Soft is a retention dressing that is particularly suitable for catheters with an injection port. It consists of white non-woven material, adhesive material and cover paper. The dressing is provided with a separate padding cushion.

Product composition:
Poliéster, polyacrylate, viscoso, polypropylene, polyethylene, glassine paper with polydimethylsiloxane coating

Intended use:
For the retention of cannulas and catheters on the skin and as a mechanical barrier to protect the puncture site.

Indications:
Not applicable.

Contraindications:
Known allergy and/or hypersensitivity to any of the product components.

Side effects:
In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions:
The insertion site must be checked regularly for erythema, irritation and inflammation.
If infection or other signs of inflammation occur (excessive warmth, oedema, erythema, pain), the decision to continue treatment must be made by healthcare professionals.
Do not use if package is damaged or accidentally opened.

Instructions for use:
Application is performed by healthcare professionals.

- Preparation:**
After thorough cleaning and disinfection of the puncture site and the surrounding skin, the dressing is applied to the dried skin.
- Applying the retention dressing (after puncture):**
Before opening the sterile barrier, check for obvious damage.
 - Open the packaging at the point indicated and remove the product in a sterile manner (Fig. 1).
 - The separately provided padding cushion can be positioned, as padding, under the cannula shaft (Fig. 2).
 - Remove the large cover paper (Fig. 3).
 - Spread the dressing out at the wings and place the round incision around the injection port. The puncture site is now covered with the wound pad (Fig. 4).
 - Now detach the 2 small cover papers from the wings and mould the dressing in place so that they slightly overlap (Fig. 5).
- Dressing change:**
The intervals at which the dressing must be changed are determined by the treating physician. These depend on the patient's individual skin condition and on the underlying medical findings.
Replace dressings which are soaked through, soiled or no longer adhere securely. Also replace the dressing if weeping pockets develop under the dressing.

Additional note:
To detach Curafix® i.v. Soft, the use of dressing-loosening fluids (based on alcohol, saturated carbohydrates or siloxanes) is recommended.

4. Disposal:
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions:
Curafix® i.v. Soft is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

In case of serious incident, contact the manufacturer and the responsible health authorities.

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