Curafix® i.v. Film

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Instructions for Use

Product description and features

Curafix® i.v. Film is a dressing for retention of cannulas and catheters, consisting of transparent film, adhesive material and cover paper. A supporting frame with holding edges mounted on the back of the film dressing serves as an application aid. The transparency of the dressing allows monitoring of the cannula insertion site and the insertion channel. The dressing can be washed with water and iodine- or octenidine-containing disinfectants based on alcohol or water. The dressing is available in 2 versions: with or without an opening for exiting cannula or catheter tubes. The dressing is available in 2 versions: with or without an opening for exiting cannula or catheter tubes.

Product composition

Polyurethane, polyacrylate, glassine paper with polydimethylsiloxane coating, polyester, kraft paper with polyethylene coating

Intended use

For the retention of cannulas and catheters on the skin and as a mechanical barrier to protect the puncture site.

Indications

Not applicable.

Contraindications

Known allergy and/or hypersensitivity to any of the product components.

Side effects

In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions

- The insertion site must be checked regularly for erythema, irritation and inflammation.
- If infection or other signs of inflammation occur (excessive warmth, oedema, erythema, pain), the decision to continue treatment must be made by healthcare professionals.
- Do not use if package is damaged or accidentally opened.

Instructions for use

Application is performed by healthcare professionals.

1. Preparation

After thorough cleaning and disinfection of the puncture site and the surrounding skin, the dressing is applied to the dried skin.

2. Use: Applying the retention dressing (after puncture)

Before opening the sterile barrier, check for obvious damage.

- a) Open the packaging at the point indicated and remove the product in a sterile manner (Fig. 1).
- b) Turn the product so that the printed side faces you. Hold the grip tabs on the two cover papers between thumb and index finger and pull the cover papers apart completely

- (Fig. 2A/2B). To make this easier, ensure that you support the film dressing a little from below with your other fingers.
- c) Turn the dressing over so that the adhesive side faces the puncture site (Fig. 3A/3B).
- d) Place the product centrally over the puncture site and press down on it with your fingers, ensuring that the cannula or catheter tubes are also fixed in place (Fig. 4A).
- e) The Curafix® i.v. Film with an opening allows exiting cannula or catheter tubes to pass through the opening provided for this purpose (Fig. 4B).
- e) Now hold the grip tabs on the supporting frame between thumb and index finger and pull the cover paper off over the applied film (Fig. 5A/5B).

3. Dressing change

The intervals at which the dressing must be changed are determined by the treating physician. These depend on the patient's individual skin condition and on the underlying medical findings. Replace dressings which are soaked through, soiled or no longer adhere securely. Also replace the dressing if weeping pockets develop under the dressing.

Additional note

To detach Curafix® i.v. Film, the use of dressing-loosening fluids (based on alcohol, saturated carbohydrates or siloxanes) is recommended.

4. Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

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