Raucocel®

- Nasentamponade
- Méchage Nasal
- Nasal Packing
- Taponamiento nasal
- Tamponamento nasal
- Tampone nasale
- Neustampon
- Næsetampon
- Nästamponad
- Nosní tamponáda
- Nosová tamponáda
- Orrtampon
- Tamponada za nos
- Рівнікó тaмпо́н
- Ninaverejooksu tampoon

Shanghai International Holding Corp. GmbH (Europe), Eiffestrasse 80, 20537 Hamburg, Germany
 Huizhou Foryou Medical Devices Co., Ltd. North Shangxia Rd. Dongjiang Hi-tech Industry Park, 516005, Huizhou, P.R.China
 Lohmann & Rauscher International GmbH & Co. KG, Westerwaldstraße 4, 56579 Rengsdorf, Germany www.Lohmann-Rauscher.com

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en Instructions for use

Product description
Raucocel® Nasal Packing is a sterile nasal packing composed of polyvinyl acetal expansive sponge and string with or without airway tube.

Properties of product
Polyvinyl acetal sponge expands rapidly and fully after absorption of water, adapts to the internal structure of nasal cavity and wound shape, provides uniform support and compression, so hemostasis can effectively be achieved.

Adverse effects
In rare instances the physicochemical conditions associated with nasal surgery, both with and without nasal packings, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include sudden fever, vomiting, diarrhea, dizziness, fainting

Product composition
Polyvinyl acetal, string, airway tube: silicone rubber.

(or near fainting when standing up) and/or a rash looking like a sunburn. Nasal packing exhibit no antimicrobial properties. Rarely, displacement and ingestion of a nasal packing may occur. The expansive sponge is non-toxic and is considered safe to pass once it is completely in the alimentary tract. The patient should be monitored in the normal manner for the ingestion of any non-toxic material. If aspiration should occur, prompt medical intervention is imperative with close patient monitoring and airway support pending the removal of the material.

Warnings and precautions
Packaging must be checked before usage. Do not use if the packaging is damaged. Keep it dry once opened. The product should not be moistened with water or water-bearing liquid before use. When the product is used for ENT department, the string must be fixed to the wing of the nose or the cheek to prevent being aspirated into the respiratory airway. Do not keep on the wound for longer than 72 hours. This product is not reusable.

Remark
In case of serious incident, contact the manufacturer and the responsible health authorities.

Sterilization
The products are sterilized by irradiation, single use only.

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