Suprasorb® CNP endo

Gebrauchsanweisung
Mode d'emploi
Instructions for use
Istruzioni per l'uso

endo

Inhalt

Deutsch	5 – 42
Français	43 – 80
English	81 – 118
Italiano	119 – 156

 Suprasorb® endo

Suprasorb® CNP endo

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Suprasorb® CNP endo

Copyright © 2023 ASSKEA GmbH, Gebesee.

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Suprasorb® CNP endo

Copyright © 2023 ASSKEA GmbH, Gebesee.

The safety of the **Suprasorb® CNP endo** therapy unit corresponds to the recognised rules of technology and the regulations of the **German Medical Devices Act (MPG)**.

The **Suprasorb® CNP endo** therapy unit bears the **CE marking CE0494**.

The **Suprasorb® CNP endo** therapy unit has been tested in accordance with the IEC 62353 standard.

The **quality management system** in place at ASSKEA GmbH has been certified in accordance with the pertinent international standards.

The **Suprasorb® CNP endo** therapy unit is a Class IIa medical suction device.

Subject to change without notice. Errors and omissions excepted.

Suprasorb® CNP endo

CONTENTS

1	Note regarding damage on delivery and customer feedback	86
2	User instructions	86
2.1	How to use these instructions	86
2.2	Symbols	87
2.2.1	Device, packaging and accessories	87
2.2.2	Display	88
2.3	Formatting conventions	88
2.4	Glossary	88
2.5	Intended purpose	90
2.5.1	Intended use	90
2.5.2	Essential functions	90
2.5.3	Applied parts	90
2.5.4	Indications	90
2.5.5	Contraindications	90
2.5.6	Restrictions on use	90
2.6	Essential safety information – IMPORTANT!	91
2.7	User requirements	93
2.8	Product liability information	94
3	Product description	94
3.1	Illustration of the entire ME system	94
3.1.1	Suprasorb® CNP endo therapy unit	94
3.2	Scope of supply	95
3.2.1	Scope of supply of the Suprasorb® CNP endo therapy unit (sale/hire)	95
3.3	Product characteristics	95
3.3.1	Disposable exudate chamber for the Suprasorb® CNP endo therapy unit	96
3.3.2	Important information about the filter system for the Suprasorb® CNP endo therapy unit	96
3.3.3	Important information on the activated charcoal filter of the Suprasorb® CNP endo therapy unit	97
3.3.4	Important information about the battery	97
3.3.5	Pressure settings	98
3.4	Warranty	98
4	Operation	99
4.1	Set-up and commissioning	100
4.1.1	Connecting the Suprasorb® CNP endo therapy unit	100
4.2	Commissioning	100
4.2.1	Positioning the Suprasorb® CNP endo therapy unit	101

4.2.2	Connecting the disposable exudate chamber to the Suprasorb® CNP endo therapy unit	102
4.2.3	Connecting a drain	102
4.3	Operating the Suprasorb® CNP endo therapy unit	103
4.3.1	Operating the unit for the first time	103
4.3.2	Operating mode	103
4.3.3	Selecting a language	104
4.3.4	Alarm delay during start-up	105
4.4	Replacing the chamber	105
4.4.1	Replacing the disposable exudate chamber	105
4.5	Decommissioning	106
5	Servicing	107
5.1	Cleaning and disinfection	107
5.1.1	Warnings	107
5.1.2	Surface of the device	107
5.1.3	Disposing of the disposable exudate chamber	108
5.2	Reusing the therapy unit	108
5.3	Servicing and maintenance	109
6	Error resolution	110
6.1	Function test	110
6.2	Troubleshooting	110
6.3	Error messages	111
7	Transport, storage and disposal	113
7.1	Transport/returns	113
7.2	Storage	113
7.3	Disposal	113
8	Technical data	114
9	EMC information	115
9.1	Electromagnetic environment in which the Suprasorb® CNP endo therapy unit may be operated	116
9.2	Dealing with electromagnetic interactions	117
9.3	Overview of all cables and transformers that can be replaced by the user	118
10	Ordering information	118
10.1	Suprasorb® CNP endo	118
11	Legal information	118

1 Note regarding damage on delivery and customer feedback

IMPORTANT:

Check for damage on delivery and remove any packaging. Report any damage within 3 days, as otherwise any claims for damage cannot be accepted.

The products described in these instructions for use are subject to continuous development and improvement. We therefore welcome any customer feedback, comments and suggestions regarding our products and their instructions for use that help improve the product, service or documentation.

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2 User instructions

2.1 How to use these instructions

Read the instructions for use through completely before you use the **Suprasorb® CNP endo** therapy unit for the first time. For more detailed information about the components and their use, please refer to the instructions for use of the respective L&R **Suprasorb® CNP endo** Kit. If you require assistance with set-up, operation and maintenance, please contact L&R (see Section 11). Report any incidents also to L&R.

Read the safety information carefully (see Section 2.6) to avoid hazardous situations.



These instructions for use are a component of the **Suprasorb® CNP endo** therapy system. As such, you should store these instructions for use somewhere where they are easily accessible at all times. If you pass the **Suprasorb® CNP endo** therapy unit on to a third party, always include the instructions for use.

2.2 Symbols



2.2.1 Device, packaging and accessories

Symbol	Meaning
	Caution indicates possible physical injury, health risk or possible material damage.
	Keep away from moisture
	Protection class II
	Humidity limit
	Air pressure limit
	Follow the instructions for use
	Follow the instructions for use
	Degree of protection: Type BF (Body Floating)
	Temperature limit
	This device must not be disposed of as domestic waste.

Symbol	Meaning
	PLEASE NOTE Note with useful information and tips.
	IP code: IP33 (see Section 2.4)
	Order number
	Serial number
	Lot number
	Date of manufacture
	Manufacturer
	Do not use if package is damaged!
	Not for reuse
	Power supply

2.2.2 Display

Symbol	Meaning
--------	---------

	Battery full
	Battery empty
	Up
	OK (on, enter)
	Key lock (symbol on display) Automatically activated during operation and can be unlocked by pressing the up and down keys at the same time.

Symbol	Meaning
--------	---------

	Battery low
	Power supply is connected.
	Down
	Cancel (off, back)
	Filter runtime reached; internal filter is required to be replaced by L&R!

2.3 Formatting conventions

Symbol	Meaning
--------	---------

•	Bullet points
1. 2.	Perform the actions in the described order.

2.4 Glossary

C

Contamination	Contamination means that bacteria and viruses from the wound exudate come into contact with the interior of the therapy unit.
---------------	---

D

Degree of protection	The degree of protection specifies the level of protection applied parts provide against electric shock. Type BF applied parts must be set up isolated from the ground and are not suitable for direct cardiac application.
----------------------	---

I

incl.	Abbreviation for including
IP33	Ingress protection/IP code The ingress protection specifies the degree of protection of the therapy unit against contact and the ingress of liquids. The Suprasorb® CNP endo therapy unit is protected against the ingress of solid foreign bodies measuring ≥ 2.5 mm in diameter and against splashing water from any direction, even when the device is positioned at an angle of up to 60°.

L

L&R	L&R is the short form used in these instructions for use for the company Lohmann & Rauscher International GmbH & Co. KG based in Rengsdorf, Germany.
-----	--

M

ME system	Abbreviation for medical electrical system
MRI	Abbreviation for magnetic resonance imaging This uses a strong magnetic field to produce imaged sections of the inside of the body that allow the organs to be assessed.

O

Oversuction	Oversuction is when wound exudate is suctioned into the interior of the therapy unit.
-------------	---

S

Servicing	Servicing at L&R includes not only cleaning and disinfection and the associated processes but also the testing and restoration of technical-functional safety.
-----------	--

W

Wound exudate	Wound exudate refers to all accumulating fluids and particles which can develop or be present in a wound. The wound exudate is suctioned out of the wound with the Suprasorb® CNP endo therapy unit and collected in a disposable exudate chamber.
---------------	---

2.5 Intended purpose

2.5.1 Intended use

The **Suprasorb® CNP endo** therapy unit is used in combination with **Suprasorb® CNP endo** Kits to create, maintain and control negative pressure and to drain body fluids in the upper and lower gastrointestinal tract.

It is typically used in inpatient settings (professional healthcare facility environment).

2.5.2 Essential functions

- generation of negative pressure (low vacuum)
- generation of volumetric flow (low flow)

2.5.3 Applied parts

The exudate chamber tube is a type BF applied part.

2.5.4 Indications

Leaks and defects in the upper and lower gastrointestinal tract:

- perforations (iatrogenic perforations, spontaneous perforations)
- anastomoses/insufficiencies

In addition, the indications for the respective **Suprasorb® CNP endo** Kit must also be adhered to.

2.5.5 Contraindications

- defects with involvement of the bronchial system (bronchus/trachea/lungs)
- direct contact of the drain with oesophageal varices
- direct contact of the drain with exposed intestinal serosa
- pre-existing coagulation disorders (risk of bleeding)
- direct contact of the drain with large blood vessels
- irreversible ischaemia in the area being treated

In addition, the contraindications for the respective **Suprasorb® CNP endo** Kit must be adhered to. Furthermore, the side effects of the respective **Suprasorb® CNP endo** Kit must also be adhered to.

2.5.6 Restrictions on use

- in medical facilities where potential equalisation is required (e.g. cardiac surgery)
- in potentially explosive atmospheres/in the MRI environment
- outdoors
- home-care setting

2.6 Essential safety information – IMPORTANT!

In addition to the safety information provided here, the specific information given for the respective **Suprasorb® CNP endo** Kits must also be adhered to.

Warning against damage caused by incorrect voltage supply

Improper use causes overvoltages to develop in the therapy unit, which can then pass to the user.

- Ensure that the power supply is set to a voltage of 100 V to 240 V AC at a power frequency of 50-60 Hz before turning on for the first time.
- Ensure that the power supply is set to a voltage of 120 V AC before turning on for the first time in UL-listed markets such as the US and Canada.
- Only operate the therapy unit with the power supply provided (see Section 8).

Health risks associated with infectious or pathogenic microbes

Infectious and pathogenic microbes in the wound exudate can pose health risks.

- Ensure that wound care is performed carefully.
- Please follow the hygiene, cleaning and decontamination instructions.

Risk of strangulation

- People can strangle themselves with the tubes or power cable, especially if the tubes or cable are unusually long.
- Ensure that no unauthorised/uninvolved people are in the vicinity of the therapy unit during suction.
- Store the therapy unit and accessories in the case provided until needed again.

Warning against damage caused by electromagnetic phenomena

Medical electrical devices are subject to special safety precautions in terms of electromagnetic compatibility and need to be installed and commissioned in accordance with the EMC information contained in the accompanying documentation! (see Section 9)

Risk of injury through improper use

- Only use the therapy unit in accordance with the intended purpose.
- Never use the therapy unit for suction purposes in the area of the trachea.
- Never use the therapy unit for chest drainage.
- When using the power supply, first ensure that the device connector on the power supply is connected to the therapy unit and only then connect the power supply to the power connection (100 V to 240 V AC).
- The power supply must be disconnected from the power connection in exactly the opposite order – remove the power supply from the power connection (100 V to 240 V AC) first and then disconnect the device connector from the therapy unit.
- Never simultaneously touch the patient and any parts of non-ME devices near the patient.

Safety defects due to incorrect accessories and spare parts

The use of accessories and spare parts not recommended by L&R can have a detrimental effect on the safety and functioning of the therapy unit. Damage caused by the use of non-recommended accessories and spare parts or improper use is not covered by any warranties. Use only original accessories and spare parts.

Warning against safety defects due to non-permissible connections of the ME system

Connecting the ME system to other equipment or devices/pieces of equipment not recommended by L&R and not specified in the instructions for use can adversely affect the safety and functioning of the therapy unit. Damage caused by the connection of non-recommended equipment or devices/pieces of equipment to the **Suprasorb® CNP endo** therapy unit or due to improper use is not covered by any warranties.

Only connect recommended original parts with the **Suprasorb® CNP endo** therapy unit.

Damage to device due to overheating

- Never cover the power supply.
- Keep the therapy unit, the power cable and the power supply away from other sources of heat.

Damage to device through improper use

- Never suction flammable, corrosive or explosive liquids or gases.
- Do not drop the therapy unit.
- Check the housing for any damage prior to use and do not operate the therapy unit if there is visible damage to the housing.
- When using the **Suprasorb® CNP endo** therapy unit, always ensure there is sufficient light to see all markings clearly.
- To ensure the proper functioning of the device, check prior to use whether all parts susceptible to wear and tear and damage are in good working order. Should this not be the case, promptly replace the parts in question.

Damage to device due to ingress of liquids

The **Suprasorb® CNP endo** therapy unit has the IP code IP33 with respect to the ingress of liquids. Nevertheless, you should protect the therapy unit from moisture.

- Do not use the therapy unit near splashing water.
- Do not use the therapy unit in rooms with high humidity, when bathing, or when showering.
- Make sure the power supply, plug and display film do not become wet.
- Never immerse the therapy unit in water or any other liquids (even when switched off).

Checking the internal power supply

As the internal battery is not automatically kept ready for use, periodically check the battery status and, if necessary, have the battery replaced by service personnel. The battery may be replaced only by authorised service personnel, as replacement by inadequately trained individuals may result in a hazard (such as overheating, fire or explosion)!

Beware of possible physiological effects and non-obvious risks

- Select the suction level depending on the patient and the medical indication to avoid any personal injury.
- Always position the therapy unit upright on a stable and smooth level surface. Ensure that the therapy unit is not knocked over or can fall down to avoid anybody being hit by the falling therapy unit.
- Other equipment, examinations and treatments may suffer interference from the therapy unit. For this reason, always pay special attention to other equipment and to examinations and treatments running in parallel to detect any interference as rapidly as possible.
- When using the therapy unit, always ensure there is sufficient light to see all markings clearly.
- Small, detached parts may be inhaled or swallowed. Therefore ensure that no unauthorised people or children are in the vicinity of the therapy unit.
- Although the materials used have been tested for their tolerability, in exceptional cases allergic reactions to accessible materials on the therapy unit and its accessories may occur. This applies particularly to contact injuries from prolonged exposure. In this case, contact a doctor immediately.

Contact of the interior of the device with liquids and solids

If liquids or solids come into direct contact with the interior of the device, the unit must be checked by L&R.

2.7 User requirements

Therapy with the **Suprasorb® CNP endo** therapy unit may only be performed by doctors who are experienced in interventional flexible endoscopy of the gastrointestinal tract. Familiarise yourself with the operation of the **Suprasorb® CNP endo** therapy unit prior to using it for the first time.

Training on the use of the **Suprasorb® CNP endo** therapy unit is provided by L&R. The product training includes, among other things, explanations of how to assemble the device and the way it functions, the use of the therapy unit, the alarm system, its cleaning and disinfection, what to do when switching patients and how to dispose of the products.

The training should be repeated regularly every 24 months.

Every participant receives a certificate as proof of participation in the training sessions.

2.8 Product liability information

The operator assumes liability for the device's functioning if:

- the therapy unit is used other than for the specified intended purpose,
- the therapy unit is not used in accordance with the instructions for use,
- the therapy unit is opened by unauthorised persons,
- assembly, set-up, add-ons, repairs or servicing are performed by unauthorised persons,
- original accessories and spare parts are not used.

Advice to the responsible organisation:

The assembly of ME systems and modifications undertaken during their expected service life require evaluation of compliance with the requirements of the applicable standards.

3 Product description

3.1 Illustration of the entire ME system

3.1.1 Suprasorb® CNP endo therapy unit


Fig. 1 Suprasorb® CNP endo

A Disposable exudate chamber with integrated exudate chamber tube

B Chamber lock

C Buttons  (on) and  (off)

D Display

E Plus  and minus  arrow buttons

F Suprasorb® CNP endo therapy unit

G Power supply port

3.2 Scope of supply

3.2.1 Scope of supply of the **Suprasorb® CNP endo** therapy unit (sale/hire)

- **Suprasorb® CNP endo** therapy unit
- power supply incl. country adapter
- these instructions for use
- **Suprasorb® CNP endo** carry bag (for use by each individual patient)
- dispatch note and country-specific attachments

3.3 Product characteristics

The **Suprasorb® CNP endo** therapy unit may only be used in combination with **Suprasorb® CNP endo** Kits.

The **Suprasorb® CNP endo** therapy unit is a lightweight, portable, battery-operated medical suction device for stationary use for the medical suction of wound exudate. It is used for suction in the low vacuum range and can be employed in both hospitals and practices as well as when transporting patients (e.g. in an ambulance).

The **Suprasorb® CNP endo** therapy unit has a suction power of max. 8 l/min (see Section 8 “Technical data”).

The **Suprasorb® CNP endo** therapy unit can be operated via the internal battery or supplied power supply, which can also be used to charge the battery. Over-temperature protection also prevents the battery from overheating by suspending charging when the temperature of the battery is too high (e.g. due to unfavourable ambient conditions).

The vacuum is generated by a maintenance-free membrane pump driven by an electric motor. Once turned on, the negative pressure pump generates negative pressure in the tubing system and disposable exudate chamber; the negative pressure is used to suction out the wound exudate. The wound exudate is routed away from the patient and collected in the disposable exudate chamber. The integrated overflow prevention system makes the therapy unit trigger the “System closed - Canister full” alarm when the disposable exudate chamber is full, and the pump is then stopped.

The **Suprasorb® CNP endo** therapy unit may only be operated with the disposable exudate chamber (system) provided.

The disposable exudate chamber supplied with the **Suprasorb® CNP endo** therapy unit is intended for single use.

Use of the carry bag

- The **Suprasorb® CNP endo** carry bag supplied with the **Suprasorb® CNP endo** therapy unit is intended for mobile use by each individual patient. Once used, it may not be used for anybody else and must be disposed of.

3.3.1 Disposable exudate chamber for the Suprasorb® CNP endo therapy unit

The disposable exudate chamber is composed of a chamber with a connected exudate chamber tube. There is a hydrophobic bacterial filter, an activated charcoal filter and a gelling agent integrated in the disposable exudate chamber. The internal filter prevents oversuction in the case of an operating error. Once liquid reaches this filter, it is no longer possible to perform suction and the error message “System closed - Canister full” appears on the display. The suction procedure is interrupted and the disposable exudate chamber needs to be changed. The activated charcoal filter in the disposable exudate chamber reduces dissemination of odours.

Gelling agent:

The gelling agent makes it possible to transport disposable exudate chambers full of wound exudate and dispose of them without leakage. Irrespective of the suction intervals, the exudate generally thickens completely after an average gelling period of 2 to 5 minutes (depending on the characteristics of the wound exudate).

The **disposable exudate chamber incl. exudate chamber tube** is intended for **single use only**. Replace the disposable exudate chamber in accordance with the respective applicable hygiene regulations, at the latest every week, before changing patients and if it becomes full.

3.3.2 Important information about the filter system for the Suprasorb® CNP endo therapy unit

The filter system of the **Suprasorb® CNP endo** therapy unit is composed of the external filter, which is integrated in the disposable exudate chamber, and the internal self-sealing filter built into the therapy unit.

The filter system provides effective protection against oversuction and contamination of the device interior.

Product lifetime and reuse

The internal filter of the filter system is not intended for reuse. To ensure consistent functioning, the internal filter must be replaced if it comes into **contact with the exudate (blockage)**, **once the filter lifetime has elapsed** (indicated as [!⚠] on the display) or when **servicing/repairing the unit**.

The internal filter must be replaced by L&R.

3.3.3 Important information on the activated charcoal filter of the **Suprasorb® CNP endo** therapy unit

An additional filter for the elimination of undesirable odours from the exhaust air compartment of the **Suprasorb® CNP endo** therapy unit is located in the exhaust air compartment of the therapy unit. This filter is a thin piece of nonwoven material coated with activated charcoal. The activated charcoal on the nonwoven material binds the odour particles from the exhaust air and neutralises them, effectively reducing the spread of odours.

Product lifetime and reuse

The activated charcoal filter of the **Suprasorb® CNP endo** therapy unit is not intended for reuse. To ensure consistent functioning, the activated charcoal filter must be replaced when **servicing/repairing the unit or at the latest every 2 years**.

The activated charcoal filter must be replaced by L&R.

3.3.4 Important information about the battery

Checking the internal power supply

- As the internal battery is not automatically kept ready for use, check the battery status and, if necessary, have the battery replaced by service personnel.
- The battery may only be replaced by L&R, as improper replacement might result in a hazard (such as excess temperatures, fire or explosion).

The battery's charge is shown in the display.

Before using the **Suprasorb® CNP endo** therapy unit for the first time, it is expressly recommended to charge the battery fully and repeat the procedure each time after the first few uses.

The **Suprasorb® CNP endo** therapy unit is equipped with a lithium ion battery, which, unlike standard rechargeable batteries, only exhibits very low self-discharge.

The **Suprasorb® CNP endo** therapy unit should be stored and charged taking the ambient conditions specified in the technical data into account, ideally at room temperature. Never store the therapy unit including battery with the battery empty!

If the therapy unit is not used for an extended period of time (approx. 10 months), the battery should be fully recharged.

Lithium ion batteries do not have a memory effect. They can and thus should be recharged each time following a successful initial charge. Avoid frequently charging the battery for short periods of time.

Even though the battery of the **Suprasorb® CNP endo** therapy unit includes safety precautions against deep discharging, the charging information should still be taken into account. The battery is also protected against overheating during charging. Should the battery temperature be exceeded during charging due to impermissible ambient conditions, the charging is temporarily interrupted to allow cooling. This is a safety measure that protects the battery.

The battery manufacturer states that the battery retains a residual capacity of more than 80% after 300 charging cycles.

3.3.5 Pressure settings

The pressure settings can be adjusted by the healthcare professional as soon as the **Suprasorb® CNP endo** therapy unit is switched on.

Pressure settings between -40 mmHg and -125 mmHg (in 5 mmHg intervals) are possible. Instructions on how to set the pressure and time values can be found in Section 4.2.

All settings on the therapy unit are set only by doctors or by healthcare professionals under the direction of a doctor.

Ensure that the **Suprasorb® CNP endo** therapy unit is equipped with a disposable exudate chamber before switching it on.

Perform early endoscopic monitoring in cases of suspected therapy failure.

3.4 Warranty

Modifications to the ME device are not permitted.

The **Suprasorb® CNP endo** therapy unit comes with a 2-year warranty. This period is neither extended nor renewed when warranty work is performed. The battery comes with a 6-month warranty. The warranty does not extend to wearing parts or consumables.

L&R is only responsible for consequences for safety, reliability and specified performance capacity if:

- original accessories and spare parts are used,
- repairs are performed by L&R,
- the **Suprasorb® CNP endo** therapy unit is used and operated in accordance with the instructions for use and not contrary to the intended purpose.

All warranty claims are null and void if

- the therapy unit is opened by unauthorised persons,
- the security seal is removed or broken,
- repairs are performed by unauthorised persons,
- modifications are made to the **Suprasorb® CNP endo** therapy unit, as in such cases the fundamental safety of the **Suprasorb® CNP endo** therapy unit can no longer be guaranteed and the functioning of the device may be impaired.

4 Operation

Risk of injury through improper use

- Only use the therapy unit in accordance with the intended purpose.
- Only use suitable suction equipment for suctioning.
- Read Sections 4.1 and 4.3 completely!

Risk to patient when lowering the pressure value

If the pressure value is lowered too much, it cannot be ruled out that a leak in the drain may not be detected and thus no alarm is triggered. The **Suprasorb® CNP endo** therapy unit has been developed to generate negative pressure in the drain and to suction off wound exudate.

Risk to patient from heavy bleeding

The therapy unit is not designed to detect heavy bleeding or pooling of blood in the wound. For this reason, it is not sufficient to rely solely on the alarm system of the **Suprasorb® CNP endo** therapy unit.

Risk to patient from severe bleeding

If severe bleeding is suspected, immediately perform endoscopic monitoring of the therapy under surgical conditions and in readiness to perform surgery.

Device failure as a result of sucked-in wound exudate

- Ensure that the disposable exudate chamber for the **Suprasorb® CNP endo** therapy unit is replaced regularly. Once the disposable exudate chamber is full, the integrated overflow prevention system is activated and the alarm “System closed - Canister full” is triggered. This suspends suction.
- Switch off the therapy unit when replacing the disposable exudate chamber.
- If the internal filter of the **Suprasorb® CNP endo** therapy unit is blocked, the therapy unit must be serviced by L&R!

Damage to device due to inadequate acclimatisation

After the devices have been exposed to ambient temperatures outside the technical specifications (see Section 8) during transport/storage, they must first become acclimatised at room temperature (approx. 20°C) for approx. 2 hours to allow use as intended.

Regularly monitoring patients can reduce the risk of deterioration of the patient's condition.

4.1 Set-up and commissioning

The following sections describe the connection, commissioning and operation of the **Suprasorb® CNP endo** therapy unit.

4.1.1 Connecting the **Suprasorb® CNP endo** therapy unit

Check the power supply incl. the power cable for any damage and replace immediately in case of damage.

Use the plug of the **Suprasorb® CNP endo** therapy unit (Section 3.1.1, Fig. 1 [G]) to connect the therapy unit to the mains using the power supply provided (type: ATM024T-W120V) for charging or operation as necessary. Ensure that you set up the therapy unit in such a way that it can be easily disconnected later.

Important information on the permissible ambient conditions during operation can be found in Section 8, “Technical data”.

Only use the supplied power supply. Start by connecting the power supply to the power socket on the **Suprasorb® CNP endo** therapy unit and then to the mains power.

4.2 Commissioning

- Read the instructions for use through completely before you use the **Suprasorb® CNP endo** therapy unit for the first time.
- Remove the therapy unit and the accessories from the packaging.
- Always place the therapy unit on a stable, level surface ensuring that it is positioned correctly (see Section 4.2.1).
- All the tubes and the power supply need to be checked for damage each time before you use the **Suprasorb® CNP endo** therapy unit. Avoid kinking the tubes when connecting the tubes. Ensure that the disposable exudate chamber and tubing are connected correctly before switching on the unit.
- Charge the battery completely before using it for the first time.
- Perform a function test (see Section 6.1)

4.2.1 Positioning the **Suprasorb® CNP endo** therapy unit

The **Suprasorb® CNP endo** therapy unit can be positioned next to the patient's bed or affixed using the **Suprasorb® CNP endo** carry bag. The carry bag is also suitable for mobile use. However, the doctor is responsible for deciding whether the patient's condition permits mobile use. In addition, the **Suprasorb® CNP endo** therapy unit can also be used on its side:



Fig. 2
*Suprasorb® CNP endo
therapy unit on its side
with the chamber connected*

Position the **Suprasorb® CNP endo** therapy unit below the wound to be treated for optimal removal of the wound exudate. While doing so, ensure that the exudate chamber tube does not hang down and is at least level with the patient.

4.2.2 Connecting the disposable exudate chamber to the **Suprasorb® CNP endo** therapy unit

- A** Disposable exudate chamber incl. exudate chamber tube
- B** Chamber lock
- C** Suction port
- D** **Suprasorb® CNP endo**
- E** Guides

Fig. 3 Connecting the disposable exudate chamber

1. Remove the disposable exudate chamber (Fig. 3[A]) from the packaging.
2. Push the chamber onto the guides (Fig. 3[E]) on the **Suprasorb® CNP endo** therapy unit until the disposable exudate chamber can be heard to click into position audibly (Fig. 3[B]).

4.2.3 Connecting a drain

Connect the exudate chamber tube of the disposable exudate chamber to the drain.

4.3 Operating the Suprasorb® CNP endo therapy unit

4.3.1 Operating the unit for the first time

1. Press the  button for 1-2 seconds to switch the **Suprasorb® CNP endo** therapy unit on. The following start screen will appear for 5 seconds:

2. The *Setup - Language* menu will appear:
3. Select the desired language using the plus and minus buttons:  

4. Press the  button to confirm your selection.
5. The overview screen will appear.

6. Perform the desired treatment (see Section 4.3.2).
7. Switch off the **Suprasorb® CNP endo** therapy unit by pressing the  button for 3 seconds.

4.3.2 Operating mode

In the continuous operation therapy mode the **Suprasorb® CNP endo** therapy unit maintains a constant negative pressure on the drain.

The value for the negative pressure can be set between -40 mmHg and -125 mmHg in 5 mmHg intervals. The default value is -125 mmHg.

1. Press the  button for 1-2 seconds to switch the **Suprasorb® CNP endo** therapy unit on. The following start screen will appear for 5 seconds:

2. The overview screen will appear.

3. You can set the prescribed negative pressure using the plus and minus buttons   and start the therapy using the  button. The therapy screen will appear. (Default target value: -125 mmHg)


4. Press the  button to stop the therapy.

- You will now be returned to the overview screen.

- Press the  button to open the valve. This allows the pressure to be actively released.

- Switch off the **Suprasorb® CNP endo** therapy unit by pressing the  button for 3 seconds.

4.3.3 Selecting a language

There is the option to select a language when starting the **Suprasorb® CNP endo** therapy unit for the first time. The language selection is saved and automatically loaded whenever the unit is started.


Proceed as follows to change the language:

- Press the  button for 1-2 seconds to switch the **Suprasorb® CNP endo** therapy unit on. The following start screen will appear for 5 seconds:

- Simultaneously press the plus and minus buttons   on the start screen. The *Setup - Language* menu will appear.

- Select the desired language using the plus and minus buttons   and confirm your selection with the  button.

4.3.4 Alarm delay during start-up

When the **Suprasorb® CNP endo** therapy unit is started, the alarm for a leaky drain “System leaky” and the alarm “System open” are delayed for 120 seconds to ensure the drain can be placed without constantly triggering any alarms.

Risk of injury due to incorrect setting

Before using the **Suprasorb® CNP endo** therapy unit on a new patient, check whether the pressure setting that has been set is suitable for this patient.

4.4 Replacing the chamber

4.4.1 Replacing the disposable exudate chamber

Health risks associated with infectious or pathogenic microbes

Infectious and pathogenic microbes in the exudate can pose health risks.

- Wear appropriate disposable gloves when replacing the disposable exudate chamber.
- Use a disposable exudate chamber on one patient only.
- Replace the disposable exudate chamber in accordance with the respective applicable hygiene regulations, at the latest every week, before every change of patient and when it becomes full.
- Any parts which come into contact with the exudate must be cleaned and disinfected or disposed of after every suction procedure.
- The exudate and the parts contaminated with it must be disposed of properly.


Fig. 4 Replacing the chamber

1. Close the drainage tube with a tube clamp or artery clamp to ensure that the negative pressure in the wound area is maintained.
2. Close the exudate chamber tube with a tube clamp (Fig. 4[1-1]).
3. Switch off the **Suprasorb® CNP endo** therapy unit or pause its operation.
4. Detach the exudate chamber tube from the drainage tube (Fig. 4[1-2]).
5. Press the lock on the top of the disposable exudate chamber, keep it depressed and remove the disposable exudate chamber from the therapy unit horizontally (Fig. 4[2]).

6. Dispose of the disposable exudate chamber and the integrated exudate chamber tube in accordance with the applicable regulations (Fig. 4[3], see Section 7.3).
7. Attach a new disposable exudate chamber to the therapy unit. Ensure that the disposable exudate chamber is securely attached to the therapy unit (Fig. 4[4]).
8. Connect the exudate chamber tube to the drainage tube.
9. Switch on the **Suprasorb® CNP endo** therapy unit or start the therapy.
10. Release the clamps on the drainage tube.

4.5 Decommissioning

Health risks associated with infectious or pathogenic microbes

Infectious and pathogenic microbes in the exudate can pose health risks.

- Wear appropriate disposable gloves when decommissioning.
- Any parts which come into contact with the exudate must be cleaned and disinfected or disposed of after every suction procedure.
- The exudate and the parts contaminated with it must be disposed of properly.

1. Switch off the therapy unit by pressing the  button for 3 seconds.
2. Disconnect the power supply from the power supply (100 V to 240 V AC) and then disconnect the device connector from the **Suprasorb® CNP endo** therapy unit.
3. For the **Suprasorb® CNP endo** therapy unit, remove the disposable exudate chamber as described in Section 4.4.1 up to and including item 5.
4. Clean the device surface in accordance with Section 5.1.2.
5. Store the therapy unit in the case provided until needed again.

5 Servicing

5.1 Cleaning and disinfection

Keep the **Suprasorb® CNP endo** therapy unit in a generally clean state.

5.1.1 Warnings

HARM TO HEALTH WHEN DEALING WITH INFECTIOUS OR PATHOGENIC MICROORGANISMS

Infectious and pathogenic microbes in the exudate can pose health risks.

- Wear appropriate disposable gloves to avoid contact with the exudate. Additional protective measures may be necessary depending on the patient's illness. Bear in mind the operator's specifications.
- Carefully perform all the steps below to achieve the required bioburden reduction.

HARM TO HEALTH WHEN DEALING WITH DISINFECTANTS

- Please observe the disinfectant manufacturer's specifications.

DEVICE DAMAGE DUE TO UNSUITABLE CLEANSERS

Cleaners and disinfectants may damage the surface of the device.

- Do not use aldehyde or acetone to clean and disinfect the device.
- With frequent cleaning and disinfection, the plastic parts of the housing may slightly change in colour. This does not have any effect on the functioning of the device.

POSSIBLE BODILY INJURY FROM ELECTRIC SHOCK

- Switch the device off before cleaning.
- Disconnect the power supply from the voltage supply and then the power supply from the device.

5.1.2 Surface of the device

Regularly clean and daily disinfect the surfaces of the therapy unit and the connected exudate chamber.

If liquids or solids come into direct contact with the interior of the therapy unit, the therapy unit must be checked by L&R.

LIMITATIONS	Stop the device when there is visible damage or when exudate enters the device. There are no other restrictions.
PREPARATION BEFORE CLEANING	Follow the instructions in Section 4.5 (steps 1-3). The device should be disconnected from the mains and placed in front of you.
CLEANING	Remove any visible soiling with a damp, but not wet, lint-free cloth. Dampen the cloth with clean water (<40°C). Applying gentle pressure, wipe the entire surface of the device.

DISINFECTION	Now disinfect the surfaces that have been cleaned. Please observe the disinfectant manufacturer's specifications. Spray all over the surface of the device. Wipe down the entire surface. Use an alcohol-based disinfectant to disinfect the device.								
	The manufacturer recommends the following disinfectants:								
	<table border="1"> <thead> <tr> <th>Surface disinfectant</th> <th>Active ingredients class</th> <th>Type</th> </tr> </thead> <tbody> <tr> <td>Incidin® Liquid (Ecolab)</td> <td>Alcohol</td> <td>ready-to-use solution</td> </tr> <tr> <td>Bacillo® AF (Bode Chemie)</td> <td>Alcohol</td> <td>ready-to-use solution</td> </tr> </tbody> </table>	Surface disinfectant	Active ingredients class	Type	Incidin® Liquid (Ecolab)	Alcohol	ready-to-use solution	Bacillo® AF (Bode Chemie)	Alcohol
Surface disinfectant	Active ingredients class	Type							
Incidin® Liquid (Ecolab)	Alcohol	ready-to-use solution							
Bacillo® AF (Bode Chemie)	Alcohol	ready-to-use solution							
DRYING	Allow the disinfected surface to dry completely.								
MAINTENANCE, MONITORING AND INSPECTION	The device should now be free of visible soiling. You can prepare the device for use again in accordance with Section 4.2.								
STORAGE	If the device is not used for a prolonged period, store in accordance with the instructions in Section 7.2.								

5.1.3 Disposing of the disposable exudate chamber

- Wear appropriate disposable gloves when replacing the disposable exudate chamber.
- Use a disposable exudate chamber on one patient only.
- Replace the disposable exudate chamber, including the exudate chamber tube, in accordance with the applicable hygiene regulations, at the latest every week, before every change of patient and when it becomes full.
- All components which come into contact with the wound exudate must be cleaned and disinfected or disposed of after every suction procedure.
- The wound exudate and all components contaminated with it must be disposed of properly.

1. Close the exudate chamber tube of the disposable exudate chamber with the tube clamp on the exudate chamber tube.
2. Dispose of the disposable exudate chamber including the exudate chamber tube in accordance with the applicable regulations. It is intended for single use only.

5.2 Reusing the therapy unit

The **Suprasorb® CNP endo** therapy unit is suitable for reuse. However, it must be cleaned and disinfected as described above before it is used on other patients.

5.3 Servicing and maintenance

The **Suprasorb® CNP endo** therapy unit must undergo an inspection by L&R every 24 months at the latest. L&R offers its customers and partners rapid and professional servicing and the performance of all necessary checks. Servicing at L&R includes not only cleaning and disinfection and the associated processes but also the testing and restoration of technical-functional and electrical safety.

WARNING: Modifying the ME device is not permitted.

Beware of any loss of performance due to ageing phenomena

Should you notice any loss of performance of the therapy unit or any other problems, please follow the troubleshooting information in Section 6 and contact L&R where applicable.

The therapy unit may only be opened and repaired by L&R.

The internal bacterial filter must be replaced if it comes into contact with the exudate (indicated as “System closed - Canister full” on the display), once the filter runtime has elapsed (indicated as [!]) on the display) or when servicing/repairing the unit.

The charcoal filter in the device must be replaced when servicing/repairing the unit, when replacing the internal bacterial filter or every 24 months at the latest.

Perform a visual inspection and function test before every use (see Section 6.1). When doing so, do not forget to also include the therapy unit accessories.

The therapy unit can be returned for repair directly to L&R.

All accessories must be cleaned and disinfected prior to their being returned. The therapy unit itself must be treated with a surface disinfectant. Please refer to the decontamination guidelines prior to dispatch (Section 7.1).

To avoid delays, notify L&R in advance of any return and affix the “Used medical device” label to the shipping box of the therapy unit.

L&R neither guarantees the error-free functioning of the **Suprasorb® CNP endo** therapy unit nor accepts liability for any material damage or bodily injury sustained if

- original L&R accessories and spare parts are not used
- the information in these instructions for use is not heeded
- assembly, set-up, modifications, add-ons and repairs are not performed by L&R
- the security seal is removed or broken.

6 Error resolution

6.1 Function test

Perform a function test without a chamber connected before you use the therapy unit for a therapy. To do so, proceed as follows:

1. Switch the therapy unit on as described in Section 4.3.
2. Start the therapy and allow the therapy unit to idle. The “System open” alarm must be triggered after 120 seconds at the latest. If the “System closed - Canister full” alarm is triggered instead, the internal filter of the **Suprasorb® CNP endo** therapy unit is blocked and needs to be replaced by L&R.
3. Acknowledge the alarm with OK.
4. Then close off the tube nozzle with one finger and restart the therapy. The “System closed - Canister full” alarm must be triggered after 5 seconds at the latest. If the alarm fails to appear even after multiple attempts, have the therapy unit checked by L&R.

6.2 Troubleshooting

Error	Possible causes	Remedy
Therapy unit does not switch on	Battery is empty.	Connect power supply.
	Display film defective.	Please contact L&R.
No suction/ no exudate flow	Exudate chamber tube closed/ blocked/kinked.	Rinse/replace exudate chamber tube. Adjust the position of the device. Check that tubes are connected correctly.
	Tube clamp closed.	Open tube clamp.
	Overflow prevention system is blocked (disposable exudate chamber).	Replace the disposable exudate chamber.
	Internal filter blocked.	Please contact L&R.
	Therapy unit is still in <i>Setup</i> mode.	Confirm selection (see Section 4.3) and start therapy.




Contact L&R if an error occurs that you cannot resolve with the measures described here.

6.3 Error messages

- The alarms are all technically triggered alarms, as they can be identified by monitoring device-specific variables.
- All the alarm messages (with the exception of “Internal error”) need to be acknowledged with OK.
- For high-priority alarm messages, the **background of the display flashes red** and the beeper sounds (3x, pause, 2x, 3x, pause, 2x) every 3 seconds.
- For low-priority alarm messages, the **background of the display turns yellow** and the beeper sounds periodically (2x) every 16 seconds.

Error message	Status/behaviour of the therapy unit	Possible cause	Remedy
	Beeper on. Pump off.	Disposable exudate chamber is not connected correctly or not at all.	Check that the connections are correct.
	Current operating mode stops.	The exudate chamber tube is not connected to the drainage tube. The drain is not in the correct position.	Check the position of the drain. Start the therapy.
	Beeper on. Current operating mode continues for 30 seconds and then stops.	The disposable exudate chamber is not connected correctly. The exudate chamber tube is not correctly connected to the drainage tube. Leak in the wound area.	Check tube connections and disposable exudate chamber. Switch off or pause the therapy unit. Inspect the wound area visually or endoscopically. Where applicable, end CO ² insufflation. Switch on the therapy unit again or start the therapy.
	Beeper on. Pump off. Current therapy mode stops.	The disposable exudate chamber is full. Exudate flow is impeded (exudate chamber tube kinked or stenosis in tube). Tube clamp closed. If the alarm is triggered even when the disposable exudate chamber is not connected, the internal bacterial filter is blocked.	Switch off or pause the therapy unit. Check the disposable exudate chamber and replace if necessary. Check tubes. Open tube clamp. Contact L&R.
	Beeper on. Pump off. Current therapy mode stops.	The disposable exudate chamber is full.	Switch off or pause the therapy unit. Check the disposable exudate chamber and replace if necessary.
		Exudate flow is impeded (exudate chamber tube kinked or stenosis in tube).	Check tubes.
		Tube clamp closed.	Open tube clamp.
		If the alarm is triggered even when the disposable exudate chamber is not connected, the internal bacterial filter is blocked.	Contact L&R.

Error message	Status/behaviour of the therapy unit	Possible cause	Remedy
	<p>Beeper on. Pump off.</p> <p>Current therapy mode stops.</p>	Battery is empty.	Connect power supply.
	<p>Beeper on. Pump off.</p> <p>Current therapy mode stops.</p>	An internal error has occurred.	Insert power supply quickly and remove again. If the error recurs 60 seconds after restart, contact L&R.
	<p>Beeper on.</p> <p>Current operating mode continues to run in background.</p>	The battery charge status is low.	Connect power supply soon.
	<p>Beeper on.</p> <p>Therapy unit stays in standby mode.</p>	<p>The therapy has not been started.</p> <p>The therapy unit has not been switched off.</p>	<p>Start the therapy.</p> <p>Switch off the therapy unit.</p>

The alarm “System leaky” can be acknowledged and muted within 30 seconds. When this is done, the **Suprasorb® CNP endo** therapy unit remains in therapy mode. If the cause of the alarm is not resolved within 60 seconds, the alarm will sound again. If the alarm is not acknowledged within the following 30 seconds, the therapy stops and the alarm will continue to be displayed. As soon as the alarm is acknowledged, the **Suprasorb® CNP endo** therapy unit will go into standby mode.

Contact L&R if an error occurs that you cannot resolve with the measures described here.

7 Transport, storage and disposal

7.1 Transport/returns

L&R offers its customers and partners rapid and professional servicing as well as the performance of all necessary checks (see Section 5).

Clean and disinfect the **Suprasorb® CNP endo** therapy unit prior to dispatch to L&R. Please refer to Section 5! Dispose of all disposable products. Include the power supply and country adapters for testing the return. Affix the supplied “Used medical device” in the shipping bag provided for this purpose inside the case!

7.2 Storage

Store the **Suprasorb® CNP endo** therapy unit in accordance with the specifications in the technical data (Section 8)!

Charge the **Suprasorb® CNP endo** therapy unit’s battery before you store the therapy unit. This ensures that it can be used at any time.

If the **Suprasorb® CNP endo** therapy unit is not used for an extended period of time (approx. 10 months), the battery should be fully recharged!

7.3 Disposal

- The device and accessories must be disposed of in accordance with the applicable regulations.
- Decontaminate the device and accessories prior to disposal.
- In accordance with EU Directive 2012/19/EC on Waste Electrical and Electronic Equipment (WEEE) and 2011/65/EC on the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (RoHS II), the device must not be disposed of as domestic waste.
- The device is registered in the German Electrical Used Equipment Register (EAR) as a small electrical device and can be taken to the local recycling centre for disposal.
- The device and accessories can also be disposed of via L&R.
- Outside the EU: Please observe any applicable national disposal regulations!

8 Technical data

Aggregate suction power*	Max. 8 l/min (low flow)	
Pressure	Max. -125 mmHg (low vacuum); Conversion factor: 1 kPa ~ 7.5 mmHg	
Chamber	Disposable exudate chamber (250 mL)	
Exudate chamber tube	PVCnoDEHP tube, 4 mm (internal dia.), length 150 cm	
Power supply	ATM024T-W120V	
Power supply rated voltage	In: AC 100 – 240 V~ / 50 – 60 Hz / 580 - 320 mA Out: DC 12 V / 2.0 A	
Maximum load current	2.0 A	
Permissible input voltage	12 V	
Power consumption at 12 V	24 W	
Degree of protection in accordance with IEC 60601-1	Type BF	
Risk class in accordance with 93/42/EEC, IX	IIa	
Protection class in accordance with IEC 60601-1	II	
IP code in accordance with IEC 60529	IP33	
CE marking	CE0494	
Sound emissions	Operation:	35 dB (A)
	High-priority alarm:	53 dB (A)
	Low-priority alarm:	52 dB (A)
Ambient conditions	Transport/storage:	-25 °C to +60 °C max. 93% humidity non-condensing
	Operation:	+5 °C to +40 °C 15% to 93% humidity non-condensing
	Air pressure:	700 hPa to 1060 hPa
Battery, rechargeable	Min. 7.2 V; lithium ion battery	
Charging time when battery empty	6-7 hours	
Energy of the battery pack	< 80 Wh	
Dimensions (H x W x D) in mm	165 x 220 x 90	
Weight (basic device)	1.2 kg	
Pressure display precision	Target pressure > -80 mmHg	max. Δ 5 %
	Target pressure ≤ -80 mmHg	max. Δ 10 %
Operating time	Continuous operation	
Runtime running on the battery	Approx. 24-48 hours, depending on motor runtime	
Expected service life	5 years	
Article number	34210	

*The information provided may vary depending on altitude, air pressure and air temperature.

9 EMC information

WARNING: The **Suprasorb® CNP endo** therapy unit may cause electromagnetic interference with equipment, examinations and treatments. For this reason, always pay special attention to other equipment and to examinations and treatments running in parallel to detect any interference as rapidly as possible.

WARNING: Portable and mobile HF communication devices (incl. peripheral devices such as antenna cable and external antennas) can interfere with medical electrical devices and should therefore not be used closer than 30 cm to any part of the **Suprasorb® CNP endo** therapy unit incl. the cable. The performance of the therapy unit may otherwise be reduced.

WARNING: The use of accessories and spare parts, transformers and cables for the **Suprasorb® CNP endo** therapy unit that are not specified or provided by L&R can cause an increase in electromagnetic emissions or a reduction in the electromagnetic immunity of the **Suprasorb® CNP endo** therapy unit and cause it to operate improperly.

Damage caused by the use of non-recommended accessories and spare parts, transformers and cables or due to improper use is not covered by any warranty.

Use only original accessories and spare parts!

WARNING: The use of the specified or supplied accessories and spare parts, transformers and cables with equipment other than the **Suprasorb® CNP endo** therapy unit can cause an increase in electromagnetic emissions or a reduction in electromagnetic immunity. Damage caused by the use of the specified or supplied accessories and spare parts, transformers and cables with other equipment or due to improper use is not covered by any warranty. Use the accessories and spare parts, transformers and cables only with the **Suprasorb® CNP endo** therapy unit!

WARNING: The use of the **Suprasorb® CNP endo** therapy unit directly beside or stacked with other equipment should be avoided, as this might result in impermissible operation. If use of the **Suprasorb® CNP endo** therapy unit directly beside or stacked with other equipment is necessary, the **Suprasorb® CNP endo** therapy unit and the other equipment should be monitored to ensure that they function correctly in the configuration used!

The **Suprasorb® CNP endo** therapy unit fulfils the requirements in accordance with IEC 60601-1-2 / EN 60601-1-2 Electromagnetic Compatibility – Medical electrical devices without any restrictions and deviations. This reduces electromagnetic interference and interactions to a minimum. Follow the information and guidelines provided to maintain the basic safety and essential functions of the **Suprasorb® CNP endo** therapy unit throughout the expected service life.

9.1 Electromagnetic environment in which the Suprasorb® CNP endo therapy unit may be operated

The **Suprasorb® CNP endo** therapy unit is designed for operation in an electromagnetic environment as described below in which the HF interference is controlled. The customer or user of the **Suprasorb® CNP endo** therapy unit should ensure that it is operated in such an environment.

The types of environment for proper operation include professional healthcare facilities and patient transport. Specific environments, such as in the vicinity of HF surgery or MRI or in which the level of EMC interference is high, are excluded.

Emission limits

Conducted and radiated HF emissions CISPR 11

Housing

Phenomenon	Test method	Immunity test level
Electrostatic discharge	IEC 61000-4-2	±8 kV contact discharge ±15 kV air discharge
Radiated HF disturbance variables	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Radiated HF disturbance variables	in accordance with frequencies and test levels in EN 60601-1-2, Table 9	see EN 60601-1-2, Table 9
Magnetic field at the supply frequency (50 Hz)	IEC 61000-4-8	30 A/m
Magnetic fields in the immediate vicinity	IEC 61000-4-39	30 kHz, 8 A/m, CW; 134.2 kHz, 65 A/m, PM: 2.1 kHz; 13.56 MHz, 7.5 A/m, PM: 50 kHz

Mains supply AC

Phenomenon	Test method	Immunity test level
Electrical fast transient disturbances/bursts	IEC 61000-4-4	±2 kV 100 kHz pulse repetition frequency
Surges	IEC 61000-4-5	±1 kV phase-phase ±2 kV phase-earth

Phenomenon	Test method	Immunity test level
Radiated HF disturbance variables	IEC 61000-4-6	3 V 0.15 MHz to 80 GHz
		6 V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage dips	IEC 61000-4-11	0 % U_T for 1/2 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°
		0 % U_T for 1 cycle 70 % U_T for 25 cycles single-phase at 0°
Voltage interruptions	IEC 61000-4-11	0 % U_T for 250 cycles

Note: U_T is the mains alternating voltage prior to application of the test level.

9.2 Dealing with electromagnetic interactions

Although electromagnetic interference and interactions with the **Suprasorb® CNP endo** therapy unit have been reduced to a minimum, electromagnetic interference between the **Suprasorb® CNP endo** therapy unit and other equipment cannot be ruled out. In any event, therefore, observe the stated requirements and the information regarding the permissible electromagnetic environment and monitor the **Suprasorb® CNP endo** therapy unit to ensure proper function and avoid adverse events for patients and users. Where applicable, select a different location for the **Suprasorb® CNP endo** therapy unit if the permissible electromagnetic environment cannot be ensured or if you have established that the **Suprasorb® CNP endo** therapy unit or adjacent equipment is not functioning properly.

As the spread of electromagnetic variables is influenced by absorption and reflection of buildings, objects and people and the field strength of stationary transmitters cannot be theoretically determined with precision beforehand, a study of the electromagnetic phenomena of the location where the **Suprasorb® CNP endo** therapy unit is to be used should be considered to determine the electromagnetic environment of the location. If unusual performance characteristics are observed, additional measures such as a modified alignment or a different location for the **Suprasorb® CNP endo** therapy unit may prove necessary. Should the essential functions of the **Suprasorb® CNP endo** therapy unit be affected by electromagnetic interference, it can be expected that volumetric flow will be reduced or absent and that either the vacuum level will be reduced or there is no vacuum build-up at all.

9.3 Overview of all cables and transformers that can be replaced by the user

Name	Specification	Maximum length
Power supply	Type: ATM024T-W120V techn. data: 100 – 240 V~ / 50-60 Hz / 580 – 320 mA (in) DC 12 V / 2.0 A (out)	4.0 m

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