

ReadyWrap®

Measuring the upper limb.

How to measure the arm:

Measure the arm in a supported and relaxed position with the elbow slightly flexed.

The measurement points refer to the following anatomical landmarks:

Circumference points:


- cA** widest width of the Palm, just below the metacarpal phalangeal joints
- cC** narrowest point of the Wrist
- cE** around the Elbow joint and into the crease of the elbow
- cG** top of the Arm, just below the axillary crease

Lengths:

- IA-C** the length between point **cA** (widest width of the palm, just below the metacarpal phalangeal joints) and point **cC** (narrowest point of the wrist)
- IC-G** the length between point **cC** (narrowest point of the wrist) and point **cG** (top of the arm, just below the axillary crease)

NB **c** = circumference, **l** = length

All length measurements should follow body contour.
All measurements to be in centimetres.



Write your measurement results in the **box** adjacent to measuring point

Delivery: to orderer directly to patient

Your reference or that of your patient's: Ms. Mr.

Client no. _____

Name/Last Name _____

Sender _____

Street _____

Phone _____

Postal Code/Place _____

Phone _____

Date of Measurement _____

Order Date _____

Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	Extra Large
ReadyWrap® Arm				
(I-C-G) (3 x lengths)	Axilla Circumference (cG) 22-29 cm Elbow Circumference (cE) 19-26 cm Wrist Circumference (cC) 13.5-17 cm	Axilla Circumference (cG) 26-36 cm Elbow Circumference (cE) 24-33 cm Wrist Circumference (cC) 15.5-20 cm	Axilla Circumference (cG) 31-43 cm Elbow Circumference (cE) 29-40 cm Wrist Circumference (cC) 17.5-23 cm	Axilla Circumference (cG) 36-50 cm Elbow Circumference (cE) 30-43 cm Wrist Circumference (cC) 19.5-26 cm
Short				
■ Beige up to 43 cm	<input type="checkbox"/> 104514	<input type="checkbox"/> 104508	<input type="checkbox"/> 104502	<input type="checkbox"/> 104520
■ Black up to 43 cm	<input type="checkbox"/> 104515	<input type="checkbox"/> 104509	<input type="checkbox"/> 104503	<input type="checkbox"/> 104521
Average				
■ Beige up to 46 cm	<input type="checkbox"/> 104516	<input type="checkbox"/> 104510	<input type="checkbox"/> 104504	<input type="checkbox"/> 104522
■ Black up to 46 cm	<input type="checkbox"/> 104517	<input type="checkbox"/> 104511	<input type="checkbox"/> 104505	<input type="checkbox"/> 104523
Tall				
■ Beige up to 49.5 cm	<input type="checkbox"/> 104518	<input type="checkbox"/> 104512	<input type="checkbox"/> 104506	<input type="checkbox"/> 104524
■ Black up to 49.5 cm	<input type="checkbox"/> 104519	<input type="checkbox"/> 104513	<input type="checkbox"/> 104507	<input type="checkbox"/> 104525

ReadyWrap® Gauntlet				
	Length (IA-C) up to 12.5 cm Wrist Circumference (cC) 13.5-17 cm Palm Circumference (cA) 18-20 cm	Length (IA-C) up to 13.5 cm Wrist Circumference (cC) 15.5-20 cm Palm Circumference (cA) 20-22 cm	Length (IA-C) up to 14.5 cm Wrist Circumference (cC) 17.5-23 cm Palm Circumference (cA) 22-24 cm	Length (IA-C) up to 15.5 cm Wrist Circumference (cC) 19.5-26 cm Palm Circumference (cA) 24-26 cm
Right ■ Beige	<input type="checkbox"/> 104498	<input type="checkbox"/> 104496	<input type="checkbox"/> 104494	<input type="checkbox"/> 104500
Right ■ Black	<input type="checkbox"/> 104499	<input type="checkbox"/> 104497	<input type="checkbox"/> 104495	<input type="checkbox"/> 104501
Left ■ Beige	<input type="checkbox"/> 104490	<input type="checkbox"/> 104488	<input type="checkbox"/> 104486	<input type="checkbox"/> 104492
Left ■ Black	<input type="checkbox"/> 104491	<input type="checkbox"/> 104489	<input type="checkbox"/> 104487	<input type="checkbox"/> 104493

Please mark (tick or cross) where applicable to indicate the garment choice and then enter the sizing, quantity and reference code.

Garment (Arm)	Length/Side	Colour	Size	Quantity	Code
Arm	<input type="checkbox"/> 43 cm	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
	<input type="checkbox"/> 46 cm				
	<input type="checkbox"/> 49,5 cm				
Gauntlet	<input type="checkbox"/> Right	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
	<input type="checkbox"/> Left				