


# ReadyWrap®

## Measuring the lower limb.

**How to measure the leg:** Ask the patient to stand, if possible. Mark and measure the leg as follows



Write your measurement results in the **box** adjacent to measuring point

### Circumference points:

- cG** 2 cm below Gluteal Fold
- cE'** 5 cm above the Knee
- cD** 2-3 cm below Popliteal Fossa
- cC** widest part of Calf
- cB** 2 cm Above Malleolus
- cY** around Instep and Heel with foot dorsi-flexed
- cA** around base of Toes

### Lengths:

- IA-Y** heel to base of Toes
- ID-B** the length between point **cB** (2 cm above the malleolus) and point **cD** (2 cm below popliteal fossa)
- IG-E'** the length between point **cE'** (5 cm above the knee) and point **cG** (2 cm below gluteal fold)

### NB

**c** = circumference, **I** = length  
 All length measurements should be longitudinal and should follow body contours.  
 All measurements to be in centimetres.

Delivery:  to orderer  directly to patient

Your reference or that of your patient's:  Ms.  Mr.

Client no. \_\_\_\_\_

Name/Last Name \_\_\_\_\_

Sender \_\_\_\_\_

Street \_\_\_\_\_

Phone \_\_\_\_\_

Postal Code/Place \_\_\_\_\_

Date of Measurement \_\_\_\_\_

Phone \_\_\_\_\_

Order Date \_\_\_\_\_

# Please tick where applicable to indicate your garment choice:

	Small	Medium	Large	Extra Large	Extra Extra Large
<b>ReadyWrap® Thigh</b>					
(I(G-E) (2 x lengths)	Groin Circumference (cG) 48-63 cm Above Knee Circumference (cE) 38-48 cm	Groin Circumference (cG) 58-73 cm Above Knee Circumference (cE) 45-58 cm	Groin Circumference (cG) 68-83 cm Above Knee Circumference (cE) 52-68 cm	Groin Circumference (cG) 78-93 cm Above Knee Circumference (cE) 59-78 cm	Groin Circumference (cG) 88-103 cm Above Knee Circumference (cE) 66-88 cm
<input type="checkbox"/> Beige up to 25 cm <input type="checkbox"/> Black up to 25 cm <input type="checkbox"/> Beige up to 30 cm <input type="checkbox"/> Black up to 30 cm	<input type="checkbox"/> 104478 <input type="checkbox"/> 104479 <input type="checkbox"/> 104480 <input type="checkbox"/> 104481	<input type="checkbox"/> 104474 <input type="checkbox"/> 104475 <input type="checkbox"/> 104476 <input type="checkbox"/> 104477	<input type="checkbox"/> 104470 <input type="checkbox"/> 104471 <input type="checkbox"/> 104472 <input type="checkbox"/> 104473	<input type="checkbox"/> 104482 <input type="checkbox"/> 104483 <input type="checkbox"/> 104484 <input type="checkbox"/> 104485	<input type="checkbox"/> 104790 <input type="checkbox"/> 104791 <input type="checkbox"/> 104792 <input type="checkbox"/> 104793
<b>ReadyWrap® Knee</b>					
(one length - 20 cm)	Above Knee Circumference (cE) 38-48 cm Below Knee Circumference (cD) 28-38 cm	Above Knee Circumference (cE) 45-58 cm Below Knee Circumference (cD) 33-43 cm	Above Knee Circumference (cE) 52-68 cm Below Knee Circumference (cD) 38-48 cm	Above Knee Circumference (cE) 59-78 cm Below Knee Circumference (cD) 43-58 cm	Above Knee Circumference (cE) 66-88 cm Below Knee Circumference (cD) 53-68 cm
<input type="checkbox"/> Beige <input type="checkbox"/> Black	<input type="checkbox"/> 104464 <input type="checkbox"/> 104465	<input type="checkbox"/> 104462 <input type="checkbox"/> 104463	<input type="checkbox"/> 104460 <input type="checkbox"/> 104461	<input type="checkbox"/> 104466 <input type="checkbox"/> 104467	<input type="checkbox"/> 104788 <input type="checkbox"/> 104789
<b>ReadyWrap® Calf</b>					
(I(D-B) (2 x lengths)	Calf Circumference (cC) 28-38 cm Ankle Circumference (cB) 18-23 cm	Calf Circumference (cC) 33-43 cm Ankle Circumference (cB) 20-28 cm	Calf Circumference (cC) 38-48 cm Ankle Circumference (cB) 25-33 cm	Calf Circumference (cC) 43-58 cm Ankle Circumference (cB) 30-38 cm	Calf Circumference (cC) 53-68 cm Ankle Circumference (cB) 35-43 cm
<input type="checkbox"/> Beige up to 30 cm <input type="checkbox"/> Black up to 30 cm <input type="checkbox"/> Beige up to 35 cm <input type="checkbox"/> Black up to 35 cm	<input type="checkbox"/> 104450 <input type="checkbox"/> 104451 <input type="checkbox"/> 104452 <input type="checkbox"/> 104453	<input type="checkbox"/> 104446 <input type="checkbox"/> 104447 <input type="checkbox"/> 104448 <input type="checkbox"/> 104449	<input type="checkbox"/> 104442 <input type="checkbox"/> 104443 <input type="checkbox"/> 104444 <input type="checkbox"/> 104445	<input type="checkbox"/> 104454 <input type="checkbox"/> 104455 <input type="checkbox"/> 104456 <input type="checkbox"/> 104457	<input type="checkbox"/> 104784 <input type="checkbox"/> 104785 <input type="checkbox"/> 104786 <input type="checkbox"/> 104787
<b>ReadyWrap® Foot</b>					
(Side)	Foot Length (I(A-Y) up to 21 cm Ankle Circumference (cB) 18-23 cm Circumference (cY) 28-33 cm	Foot Length (I(A-Y) up to 23 cm Ankle Circumference (cB) 20-28 cm Circumference (cY) 30-38 cm	Foot Length (I(A-Y) up to 26 cm Ankle Circumference (cB) 25-33 cm Circumference (cY) 35-43 cm	Foot Length (I(A-Y) up to 29 cm Ankle Circumference (cB) 30-38 cm Circumference (cY) 40-48 cm	
Right <input type="checkbox"/> Beige Right <input type="checkbox"/> Black Left <input type="checkbox"/> Beige Left <input type="checkbox"/> Black	<input type="checkbox"/> 104436 <input type="checkbox"/> 104437 <input type="checkbox"/> 104428 <input type="checkbox"/> 104429	<input type="checkbox"/> 104434 <input type="checkbox"/> 104435 <input type="checkbox"/> 104426 <input type="checkbox"/> 104427	<input type="checkbox"/> 104432 <input type="checkbox"/> 104433 <input type="checkbox"/> 104424 <input type="checkbox"/> 104425	<input type="checkbox"/> 104438 <input type="checkbox"/> 104439 <input type="checkbox"/> 104430 <input type="checkbox"/> 104431	
<b>ReadyWrap® Toes</b>					
(cA) (Side)	Circumference (cA) 18.0-20.4 cm	Circumference (cA) 20.5-23.5 cm	Circumference (cA) 23.6-26.0 cm	Circumference (cA) 26.1-28.5 cm	Circumference (cA) 28.6-31.0 cm
Right <input type="checkbox"/> Beige Right <input type="checkbox"/> Black Left <input type="checkbox"/> Beige Left <input type="checkbox"/> Black	<input type="checkbox"/> 104416 <input type="checkbox"/> 104417 <input type="checkbox"/> 104406 <input type="checkbox"/> 104407	<input type="checkbox"/> 104414 <input type="checkbox"/> 104415 <input type="checkbox"/> 104404 <input type="checkbox"/> 104405	<input type="checkbox"/> 104412 <input type="checkbox"/> 104413 <input type="checkbox"/> 104402 <input type="checkbox"/> 104403	<input type="checkbox"/> 104418 <input type="checkbox"/> 104419 <input type="checkbox"/> 104408 <input type="checkbox"/> 104409	<input type="checkbox"/> 104420 <input type="checkbox"/> 104421 <input type="checkbox"/> 104410 <input type="checkbox"/> 104411
<b>ReadyWrap® Extender Strap</b>					
	One size fits all				
<input type="checkbox"/> Beige <input type="checkbox"/> Black	<input type="checkbox"/> 104740 <input type="checkbox"/> 104741				

Please mark (tick or cross) where applicable to indicate the garment choice and then enter the sizing, quantity and reference code.

Garment (Right Leg)	Length/ Side	Colour	Size	Quantity	Code
Thigh	<input type="checkbox"/> 25 cm <input type="checkbox"/> 30 cm	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Knee		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Calf	<input type="checkbox"/> 30 cm <input type="checkbox"/> 35 cm	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Foot		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Toe		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Extender Strap		<input type="checkbox"/> Beige <input type="checkbox"/> Black			

Garment (Left Leg)	Length/ Side	Colour	Size	Quantity	Code
Thigh	<input type="checkbox"/> 25 cm <input type="checkbox"/> 30 cm	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Knee		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Calf	<input type="checkbox"/> 30 cm <input type="checkbox"/> 35 cm	<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Foot		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Toe		<input type="checkbox"/> Beige <input type="checkbox"/> Black			
Extender Strap		<input type="checkbox"/> Beige <input type="checkbox"/> Black			