



People.Health.Care.

VENOSAN® Compression stockings

Product program.



Product overview

Medical compression stockings and pantyhose.

VENOSAN® 3000 COTTON

NATURAL LOOK

- Durable
- 100% cotton on the skin
- Absorbs moisture
- Easy to don

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- After vein surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema



Composition

- CCL 2 49% cotton
27% nylon polyamide
24% elastane
- CCL 3 40% cotton
27% nylon polyamide
33% elastane

VENOSAN® 4000

CLIMATE EFFECT

- Climate effect
- Semi-transparent
- Comfortable
- Highly elastic

Indications

- CVI (chronic venous insufficiency) 1
- Varicosis during pregnancy
- Milder symptoms
- After vein surgery



Composition

- 71% TACTEL® polyamide
29% elastane

VENOSAN® 5000

THE MAGIC OF SILVER

- Odour-inhibiting
- Anti-static
- Temperature-regulation
- Anti-bacterial

Indications

- CVI (chronic venous insufficiency) 1–2
- After vein surgery
- Mild skin conditions
- Lipoedema stage 1



Composition

- 40% TACTEL® polyamide
30% elastane
20% cotton
10% polyamide/silver (X-Static®)

VENOSAN® 6000 CLASSIC

- Maximum sheerness
- Breathable
- Sits perfectly
- Ease of care

Indications

- CVI (chronic venous insufficiency) 1-2
- Varicose veins
- Varicosis during pregnancy
- After vein surgery



Composition
70% nylon polyamide
30% elastane

VENOSAN® 7000 NATURAL TOUCH

- High stiffness (resistance)
- Solide
- Skin-friendly
- Temperature regulation

Indications

- CVI (chronic venous insufficiency) 2-3
- Lipoedema stage 2-3
- After vein surgery
- After thrombosis
- PTS (post-thrombotic syndrome)
- Mild lymphoedema
- Marked orthostatic stress
















Composition
62 % TACTEL® polyamide
26 % elastane
8 % cotton
4 % Lyocell/SeaCell®


Colours and styles

Medical compression stockings and pantyhose.

Colours

Line	Colours	Actual product colours may vary slightly from the colours in this catalogue due to the printing process.
VENOSAN® 3000 COTTON NATURAL LOOK	Beige	
VENOSAN® 4000 CLIMATE EFFECT	Mexico	
	Marokko	
	Sand	
	Black	
VENOSAN® 5000 THE MAGIC OF SILVER	Mexico	
	Marokko	
	Sand	
	Black	
VENOSAN® 6000 CLASSIC	Beige	
	Black	
VENOSAN® 7000 NATURAL TOUCH	Beige	
	Black	


Overview styles



AD
Below knee
stocking




ADH
Below knee
stocking with
dotted top band




AG
Thigh-length
stocking




AGH
Thigh-length
stocking with
dotted top band




AGH
Thigh-length
stocking with
lace top band




AGG
Thigh-length
stocking with belt



AT
Pantyhose




ATU Materna
Maternity
pantyhose



ATH
Pantyhose
for men

Product lines

Medical compression stockings and pantyhose.

											
Line	CCL	Pressure at ankle (mmHg)	AD	ADH	AG	AGH	AGH	AGG	AT	ATU Materna	ATH
VENOSAN® 3000 COTTON NATURAL LOOK	CCL 2	23–32	●	●	○	●	○	○	○		○
	CCL 3	34–46	●	○	○	○	○	○	○		○
VENOSAN® 4000 CLIMATE EFFECT	CCL 1	18–21	●	●	○	●	●	○	●	●	
	CCL 2	23–32	●	●	○	●	●	●	●	●	
VENOSAN® 5000 THE MAGIC OF SILVER	CCL 1	18–21	●	●	○	●	●	●	●	○	
	CCL 2	23–32	●	●	○	●	●	●	●	○	
	CCL 3	34–46	●	○	○	○	○	○	○	○	
VENOSAN® 6000 CLASSIC	CCL 1	18–21	⊙			⊙			⊙		
	CCL 2	23–32	⊙			⊙			⊙		
VENOSAN® 7000 NATURAL TOUCH	CCL 1	18–21	⊙								
	CCL 2	23–32	●	●	○	●	○	●	○		○
	CCL 3	34–46	●	●	○	●	○	○	○		○

- ⊙ available in standard sizes
- available in standard and custom sizes
- available in custom sizes only

Product overview

Medical support stockings.

VENOSAN® LEGLINE 20

ELEGANT LUXURY

- For non-heavy legs in everyday situations
- Retains shape, ensures maximum freedom of movement
- Sheer look
- Optimal wearing comfort

Use

- To prevent varicose veins and heavy legs
- For people who sit or stand a lot at work
- Long journeys by plane, train, and car



20 mmHg
140 den

Composition
78 % polyamide
22 % LYCRA® elastane

VENOSAN® LEGLINE 30

ELEGANT LUXURY

- For non-heavy legs in everyday situations
- Retains shape, ensures maximum freedom of movement
- Sheer look
- Optimal wearing comfort

Use

- To prevent varicose veins and heavy legs
- For people who sit or stand a lot at work
- Long journeys by plane, train, and car



30 mmHg
280 den

Composition
AD / AGH:
70 % polyamide
30 % LYCRA® elastane
AT:
78 % polyamide
22 % LYCRA® elastane

Travel socks

VENOSAN® JET LEGS®

TRAVEL WITH COMFORT SOCKS

- Gradually decreasing compression (15-20 mmHg)
- Promotes blood circulation
- Efficient thrombosis prophylaxis

Use

- Wear for plane, train, and car travel
- Long journeys



Composition
86 % nylon polyamide
14 % LYCRA® elastane

VENOSAN® COTTON SUPPORT

TRAVEL & LEISURE

- For non-heavy legs for everyday situations and travel
- Promotes blood circulation
- Efficient thrombosis prophylaxis
- 100% cotton on the skin
- Can be worn with jeans and suits

Use

- For people who sit or stand a lot at work
- Long journeys by plane, train, and car



15–20
mmHg

Composition
65 % polyamide
20 % cotton
15 % elastane

VENOSAN® SILK SUPPORT

TRAVEL & LEISURE

- For non-heavy legs for everyday situations and travel
- Promotes blood circulation
- Efficient thrombosis prophylaxis
- Soft, smooth silk fibres on the skin
- Temperature-equalising

Use

- For people who sit or stand a lot at work
- Long journeys by plane, train, and car
- For sensitive and dry skin



15–20
mmHg

Composition
65 % polyamide
20 % elastane
15 % silk

Anti-embolism stockings

VENOSAN® AES

- Perfect fit thanks to the anatomical shape
- Colour coding for simple and reliable identification of the size
- Top inspection window for safety

Indications

- To prevent thrombosis and embolism before, during and after operations
- General immobility, confinement to bed
- Before, during and after childbirth




Composition
AD / AGG:
91 % nylon polyamide
9 % LYCRA® elastane
AG:
89 % nylon polyamide
11 % LYCRA® elastane



People.Health.Care.



Lohmann & Rauscher GmbH & Co. KG · Westerwaldstraße 4 · D-56579 Rengsdorf
Lohmann & Rauscher GmbH · Johann-Schorsch-Gasse 4 · A-1141 Vienna

VENOSAN® is a registered trademark of
SWISSLASTIC AG St. Gallen, Switzerland – an  Company

2002669 0220 e